

Strong Past. Promising Future.

Ask for your free
Mammotome® trial

The accurate and safe ultrasound biopsy

A Fast and Easy Ultrasound Breast Biopsy

EXperience the design

- **EX**perience the choice of 11-Gauge or 8-Gauge probe
- Comfort fit with ergonomic contours for various hand positions
- Easy access to tissue specimens
- Simple introduction of tissue marker using special ports
- Visual indication of penetration depth with needle markings

EXpedite biopsy procedure

- Quick and easy 3-step setup
- Improved control of the device with a "cable-free" handpiece
- Smooth and controlled insertion with a bladed tip
- Rapid tissue acquisition using a high-speed cutter



Ideal option for minimally invasive biopsy, accurate tissue sampling and removal of imaged evidence

Mammotome®
EX


Mammotome®
www.mammotome.com




Strong Past. Promising Future.

The accurate and safe ultrasound biopsy

Enlarged detail
of aperture
of probe tip



Enlarged detail
of needle
markings



Enlarged detail
of open tissue
collection



MammoMARK™
placement



Ergonomically designed handpiece



Mammotome® EX



DEVICOR MEDICAL PRODUCTS

Devicor Medical UK Ltd • 20-22 Bedford Row • London WC1R 4JS • United Kingdom
info.uk@mammotome.com • www.mammotome.com

Information for European countries outside of the UK by:
DEVICOR MEDICAL GERMANY GmbH • Nordport Towers • Südportal 3 • D-22848 Norderstedt, Germany
Tel.: +49 (0) 40 - 59 35 59 10 • Fax: +49 (0) 40 - 500 98 940 • info-europe@mammotome.com
www.mammotome.com

Mammotome®
www.mammotome.com